

One-Pager

Needles

Zero Residual™ Needles



REF: ZRN3309

33G | 9mm



REF: ZRN3009

30G | 9mm



REF: ZRN30G

30G | 13mm



REF: ZRN30G-AT

30G | 16mm

Zero Residual™ provides sterile single-use injection needles with low dead volume (LDV/residual volume/dead space). All needles are precision-engineered using high quality steel, to ensure optimal sharpness and a smooth injecting experience.

To achieve optimal low dead volume, it is recommended to use the Zero Residual™ needle with the Zero Residual™ syringe.

The needles are securely protected by a needle shield and individually packed in sterile blisters. This ensures optimal conditions for ophthalmic injections, maintaining the needle's cleanliness right up to the moment of use.

Compatibility

- USP 789 compliant
- CE & FDA Approved
- Each batch is tested to ensure optimal sharpness
- Compliant with injection assistant tools